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- Omuro A, DeAngelis LM. Glioblastoma and other malignant gliomas: a clinical review. JAMA. 2013; 310:1842-50.
- Kesari S. Understanding glioblastoma tumor biology: the potential to improve current diagnosis and treatments. Semin Oncol. 2011;38 (Suppl 4):S2-10.
- Mrugala MM. Advances and challenges in the treatment of glioblastoma: a clinician's perspective. Discov Med. 2013;15:221-30.
- Stupp R, Hegi ME, Mason WP, van den Bent MJ, Taphoorn MJ, Janzer RC, et al. European Organisation for Research and Treatment of Cancer Brain Tumour and Radiation Oncology Groups. National Cancer Institute of Canada Clinical Trials Group. Effects of radiotherapy with concomitant and adjuvant temozolomide

- versus radiotherapy alone on survival in glioblastoma in a randomised phase III study: 5-year analysis of the EORTC-NCIC trial. *Lancet Oncol.* 2009;10:459-66.
5. Koshy M, Villano JL, Dolecek TA, Howard A, Mahmood U, Chmura SJ, et al. Improved survival time trends for glioblastoma using the SEER 17 population-based registries. *J Neurooncol.* 2012;107:207-12.
 6. Darefsky AS, King JT Jr, Dubrow R. Adult glioblastoma multiforme survival in the temozolomide era: a population-based analysis of Surveillance, Epidemiology, and End Results registries. *Cancer.* 2012;118:2163-72.
 7. Thon N, Kreth S, Kreth FW. Personalized treatment strategies in glioblastoma: MGMT promoter methylation status. *Onco Targets Ther.* 2013;6:1363-72.
 8. Wong ET, Hess KR, Gleason MJ, Jaeckle KA, Kyritsis AP, Prados MD, et al. Outcomes and prognostic factors in recurrent glioma patients enrolled onto phase II clinical trials. *J Clin Oncol.* 1999;17:2572-8.
 9. Brem H, Piantadosi S, Burger PC, Walker M, Selker R, Vick NA, et al. Placebo-controlled trial of safety and efficacy of intraoperative controlled delivery by biodegradable polymers of chemotherapy for recurrent gliomas. The Polymer-brain Tumor Treatment Group. *Lancet.* 1995;345:1008-12.
 10. Iwamoto FM, Abrey LE, Beal K, Gutin PH, Rosenblum MK, Reuter VE, et al. Patterns of relapse and prognosis after bevacizumab failure in recurrent glioblastoma. *Neurology.* 2009;73:1200-6.
 11. Quant EC, Norden AD, Drappatz J, Muzikansky A, Doherty L, LaFrankie L, et al. Role of a second chemotherapy in recurrent malignant glioma patients who progress on bevacizumab. *Neuro Oncol.* 2009;11: 550-5.
 12. Instructions for Use. NovoTTF-100A System. March 3, 2012.
 13. Stupp R, Wong ET, Kanner AA, Steinberg D, Engelhard H, Heidecke V, et al. NovoTTF-100A versus physicians' choice chemotherapy in recurrent glioblastoma: A randomised phase III trial of a novel treatment modality. *Eur J Cancer.* 2012;48:2192 -202.
 14. Kirson ED, Dbaly V, Tovarys F, Vymazal J, Soustiel JF, Itzhaki A, et al. Alternating electric fields arrest cell proliferation in animal tumor models and human brain tumors. *PNAS.* 2007;104:10152-7.
 15. Moritz AR, Henriques FC Jr. Studies of thermal injury. *Am J Pathol.* 1947;23:695-720.
 16. NovoTTF-100A System – P100034. Available from: <http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/DeviceApprovalsandClearances/Recently-ApprovedDevices/ucm254480.htm>.
 17. Kirson ED, Gurvich Z, Schneiderman R, Dekel E, Itzhaki A, Wasserman Y, et al. Disruption of cancer cell replication by alternating electric fields. *Cancer Res.* 2004;64:3288-95.
 18. Pless M, Weinberg U. Tumor treating fields: concept, evidence and future. *Expert Opin Investig Drugs.* 2011;20:1099-106.
 19. Davies AM, Weinberg U, Palti Y. Tumor treating fields: a new frontier in cancer therapy. *Ann N Y Acad Sci.* 2013;1291:86-95.
 20. Lee SX, Tunkyi A, Wong ET, Swanson KD. Mitosis interference of cancer cells by NovoTTF-100A causes decreased cellular viability [abstract]. *Neuro-Oncology (Meeting Abstracts).* 2012;14(Suppl 6):iv9.
 21. Lee SX, Wong ET, Swanson KD. Mitotic interference of cancer cells during anaphase by electric field from NovoTTF-100A [abstract]. *Neuro-Oncology (Meeting Abstracts).* 2011;13(Suppl 3):iii13-iii14.
 22. Kirson ED, Schneiderman RS, Dbaly V, Tovarys F, Vymazal J, Itzhaki A, et al. Chemotherapeutic treatment efficacy and sensitivity are increased by adjuvant alternating electric fields (TTFields). *BMC Med Phys.* 2009;9:1.
 23. Edgerly M, Fojo T. Is there room for improvement in adverse event reporting in the era of targeted therapies? *J Natl Cancer Inst.* 2008;100:240-2.
 24. Chen AP, Setser A, Anadkat MJ, Cotliar J, Olsen EA, Garden BC, et al. Grading dermatologic adverse events of cancer treatments: the Common Terminology Criteria for Adverse Events Version 4.0. *J Am Acad Dermatol.* 2012;67:1025-39.

25. NCI CTC Version 4 .03, 2010 available at http://evs.nci.nih.gov/ftp1/CTCAE/CTCAE_4.03_2010-06-14_QuickReference_8.5x11.pdf.
26. Proksch E, Brandner JM, Jensen JM. The skin: an indispensable barrier. *Exp Dermatol*. 2008;17:1063-72.
27. Nemes Z, Steinert PM. Bricks and mortar of the epidermal barrier. *Exp Mol Med*. 1999;31:5-19.
28. Madison KC. Barrier function of the skin: "la raison d'être" of the epidermis. *J Invest Dermatol*. 2003;121:231-41.
29. Moens E, Veldhoen M. Epithelial barrier biology: good fences make good neighbours. *Immunology*. 2012;135:1-8.
30. Rheinwald JG, Green H. Epidermal growth factor and the multiplication of cultured human epidermal keratinocytes. *Nature*. 1977;265:421-4.
31. Cruickshank CN. The microanatomy of the epidermis in relation to trauma. *J Tissue Viability*. 2006;16:16-9.
32. Cowley K, Vanoosthuyze K. Insights into shaving and its impact on skin. *Br J Dermatol*. 2012;166(Suppl. 1): 6-12.
33. Wolf R, Orion E, Ruocco E, Baroni A, Ruocco V. Contact dermatitis: facts and controversies. *Clin Dermatol*. 2013;31:467-78.
34. Voegeli D. Moisture-associated skin damage: an overview for community nurses. *Br J Community Nurs*. 2013;18:6,8,10-12.
35. Vlassova N, Han A, Zenilman JM, James G, Lazarus GS. New horizons for cutaneous microbiology: the role of biofilms in dermatological disease. *Br J Dermatol*. 2011;165:751-9.
36. Hardy MA. The biology of scar formation. *Phys Ther*. 1989;69:1014-24.
37. Chen HX, Cleek JN. Adverse effects of anticancer agents that target the VEGF pathway. *Nat Rev Clin Oncol*. 2009;6:465-77.
38. Raab WP. Photodamaged skin: a medical or a cosmetic concern? *J Int Med Res*. 1990;18(Suppl 3):2C-7C.
39. Malak JA, Kibbi AG. Revised terminology in dermatology: a call for the new millennium. *Arch Dermatol*. 2001;137:93-4.
40. Mark BJ, Slavin RG. Allergic contact dermatitis. *Med Clin North Am*. 2006;90:169-85.
41. Cohen DE, Moore MM. Occupational skin disease. In: Rom WN, Markowitz SB, editors. *Environmental and occupational medicine*. 4th ed. Philadelphia, PA: Lippincott-Williams & Wilkins; 2007. p.617.
42. Rajagopalan R, Anderson R. Impact of patch testing on dermatology-specific quality of life in patients with allergic contact dermatitis. *Am J Contact Dermat*. 1997;8:215-21.
43. Nakada T, Hostynek JJ, Maibach HI. Use tests: ROAT (repeated open application test)/PUT (provocative use test): an overview. *Contact Dermat*. 2000;43:1-3.
44. Nosbaum A, Vocanson M, Rozieres A, Hennino A, Nicolas JF. Allergic and irritant contact dermatitis. *Eur J Dermatol*. 2009;19:325-32.
45. Grice EA, Segre JA. The skin microbiome. *Nat Rev Microbiol*. 2011;9:244-53.
46. Stenn KS, Paus R. Controls of hair follicle cycling. *Physiol Rev*. 2001;81:449-94.
47. Elzinga G, Chung A, Wong ET. Safety Analysis of Bevacizumab Plus NovoTTF-100A in Patients with Recurrent Malignant Gliomas. [abstract]. *Neuro-Oncology (Meeting Abstracts)*. 2013;15(Suppl 3): iii134.
48. Lyon CC, Smith AJ, Griffiths CE, Beck MH. The spectrum of skin disorders in abdominal stoma patients. *Br J Dermatol*. 2000;143:1248-60.
49. Szymanski KM, St-Cyr D, Alam T, Kassouf W. External stoma and peristomal complications following radical cystectomy and ileal conduit diversion: a systematic review. *Ostomy Wound Manage*. 2010;56:28-35.
50. Agarwal S, Ehrlich A. Stoma dermatitis: prevalent but often overlooked. *Dermatitis*. 2010;21:138-47.

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図 3. ヒトの皮膚構造及び細胞集団の概略図。皮膚は、表皮、真皮、皮下組織という3つの主要な層で構成されている。皮膚組織を構成する細胞および構造物によって、脱水に対する有効なバリアが維持され、紫外線放射 (UVR) や病原性微生物などの脅威から保護される。血管及びリンパ管により免疫細胞の皮膚への出入りが可能であるため、皮膚の細胞集団は皮膚の炎症系及び免疫系の要求に応じて、常に流動的に行動する。Macmillan Publishers Ltd.: Nature の許可を得て転載 (MacNeil S. Progress and opportunities for tissue-engineered skin. Nature. 2007;445:874-80, copyright 2007)。